
## **HEM - THE STORY OF YOUR LABORATORY SINCE 1988**



HEM LTD

MOSCOW • RUSSIA

www.hemltd.ru

## **EXPORT OFFER**

## FLEXISTRIP™ SPREADING SPATULAS



- Designed for safe, easy and error-free manual blood smearing.
- Made of special grade polystyrene, which eliminates cuts to hands and fingers.
- Two tapered working ends with microfiber allow making wide and narrow monolayer smears ideal for microscopy.
- Flexible design protects blood cells from possible destruction.
- Best for malaria screening and routine blood analysis.

Cat #	Min Qty	Price, EUR
511001		

# READY-TO-USE CULTURE MEDIA IN PETRI DISHES FOR MICROBIOLOGY DIAGNOSIS



- The modern technology of production and bottling of culture media.
- High-quality imported and domestic ingredients and raw materials
- Defibrinated sheep blood from E&O Labs
- Multi-stage QC
- Fill automation
- Clean-room production
- ISO 13485:2003 certified
- All produced culture media is registered in Russia as medical products (IVD).

## FECAL PARASITE CONCENTRATORS KITS



- Fecal Parasite Concentrator is the cleanest, most advanced device for detection of helminth eggs and larvae, protozoan cysts, coccidian oocysts, and microsporidian spores.
- Easy to use device provides sample sedimentation for microscopic examination, and requires the shortest processing time.
- The isolation method is a modern modification of the well-known manual formalin-ether concentrating method.

#### Features

- Closed two tube system (mixing and filter tubes with all the necessary solutions).
- · Simple, safe and efficient.
- Minimal odor, no mess.
- Standardization of process.

#### Kit content:

•	Fecal Concentrator mixing tube assembly	40 pcs/pkg
•	Fecal Concentrator filter tube assembly	40 pcs/pkg
•	Ethyl acetate bottle (40 ml)	1 pc/pkg
•	User manual	1 pc/pkg

Cat #	Min Qty	Price, EUR
580050		

## **ENTEROBIASIS SAMPLING & TEST KIT**

- Helminthiasis is the wide-spread problem for developing countries, and one of the common parasites is E. vermicularis, which is very difficult to diagnose by the conventional flotation method.
- The most suitable technique for diagnosis *E. vermicularis eggs* is the microscopic examination of a perianal applicator (glass stick or others).
- This kit uses the modified classic technique where a glue-coated spatula is used instead of a glass stick.
- The kit consists of a transport container with a labeling sticker and a spatula made of transparent polystyrene, coated with water-based glue. Transparent spatula is applied to the perianal area and can be used for microscopic examination directly.

Cat #	Min Qty	Price, EUR
580013		

## THERMAL BAGS | THERMALLY INSULATED CONTAINERS

- Ideal for transportation of any suitable items including biomedical samples with guaranteed thermal insulation.
- Capacity from 6 litres up to 48 litres.
- Can be used with cold packs or hot packs
- Lightweight and comfortable to carry
- The inner insulated chamber can be easily removed for cleaning.

Cat #	Min Qty	Price, EUR
590006 590012 590024 590036 590048		

## **AUTOCLAVE BAGS WITH STERILIZATION INDICATORS**



- Made of special grade polypropylene.
- Comply with all autoclaving protocols.
- Capacity from 8 litres up to 104 litres.
- Flexo printed with thermo-sensitive sterilization indicator.
- Safe and easy-to-use.
- Environmentally friendly due to elimination of incineration process.

Cat #	Min Qty	Price, EUR
534008		
534018		
534028		
534038		
534048		

### FOR MORE INFORMATION PLEASE CONTACT

Mr. Evgeniy Dimitko, Commercial Director

Tel: +7 495 612 43 12 ext. 112 E-mail: e.dimitko@hemltd.ru

**Ms. Irina Chirochkina,**Deputy Chief of Logistics Department

Tel: +7 495 612 43 12 ext. 221 E-mail: logist@hemltd.ru

Address: 1/12, UI. 8 Marta, Moscow, Russia Web: www.hemltd.ru